


Nu-Wave Oven vs. Conventional Oven Cooking Time Comparison


(10 minute preheating time to 350F included for conventional oven)
Tests Conducted from March 1 to March 30, 2005 by SILLIKER LABORATORIES



Nu-Wave Oven saves more than 50% cooking time over a conventional oven.

Nu-Wave Oven vs. Conventional Oven
Natural Juice Content Comparison in a Hamburger
(Weight 83 g before cooking)

Tests Conducted from March 1 to March 30, 2005 by SILLIKER LABORATORIES



Nu-Wave Oven vs. Conventional Oven Fat Content Comparison in a Hamburger

Tests Conducted from March 1 to March 30, 2005 by SILLIKER LABORATORIES

